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HumanKaryoMap-12 芯片，[www.illumina.com/products/by-type/clinical-research-products/human-karyomap.html](http://www.illumina.com/products/by-type/clinical-research-products/human-karyomap.html)

iScan 系统，[www.illumina.com/systems/array-scanners/iscan.html](http://www.illumina.com/systems/array-scanners/iscan.html)

BlueFuse-Multi 软件，[www.illumina.com/clinical/clinical\\_informatics/bluefuse.html](http://www.illumina.com/clinical/clinical_informatics/bluefuse.html)

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1. Zhang S, Lei C, Wu J, et al. The establishment and application of preimplantation genetic haplotyping in embryo diagnosis for reciprocal and Robertsonian translocation carriers *BMC Med Genomics*. 2017; 10(1):60. doi: 10.1186/s12920-017-0294-x.
2. Weier HU, Munn S, Fung J. Patient-specific probes for preimplantation genetic diagnosis of structural and numerical aberrations in interphase cells. *J Assist Reprod Genet*. 1999; 16:182–191.
3. Lersch RA, Fung J, Munné S, et al. Case-Specific, Breakpoint-Spanning DNA Probes for Analysis of Single Interphase Cells. *Genet Test*. 2000; 4:273–278.
4. Hu L, Cheng D, Gong F, et al. Reciprocal Translocation Carrier Diagnosis in Preimplantation Human Embryos. *EBioMedicine*. 2016; 14:139–147.
5. Alfarawati S, Fragouli E, Colls P, et al. First births after preimplantation genetic diagnosis of structural chromosome abnormalities using comparative genomic hybridization and microarray analysis. *Hum Reprod*. 2011; 26:1560–1574.
6. Zhao W, Wu M, Chen F, et al. Robertsonian Translocations: An Overview of 872 Robertsonian Translocations Identified in a Diagnostic Laboratory in China. *PLoS One*. 2015; 10(5): e0122647.

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© 2018 Illumina, Inc. 保留所有权利。所有商标均为 Illumina 公司或其各自所有者的财产。